

CaImAn - NeuroAnalysis

Federica Fabris s1128263

May 2024

1. Deconvolution of Calcium Traces

What is meant by 'deconvolution' of the calcium traces?

Deconvolution of calcium traces refers to the process of estimating the underlying neural activity (spike train) from the observed calcium fluorescence signals. Calcium imaging is an indirect measure of neural activity where the influx of calcium ions in neurons, which occurs when neurons fire action potentials, is detected via fluorescence. The recorded signal is typically a smoothed and delayed version of the actual spiking activity due to the kinetics of the calcium indicator and the imaging system.

What is the aim of this procedure?

The primary aim of deconvolution is to recover the timing and magnitude of the neural spikes from the noisy and temporally blurred calcium signals. This is crucial because direct interpretation of the raw calcium fluorescence data can be challenging due to noise, slow dynamics of calcium indicators, and overlapping signals from multiple spikes. Deconvolution helps in providing a clearer, more accurate representation of neural activity, which is essential for studying brain function and neural circuits.

2. Dealing with Motion Artifacts in CaImAn

How does CaImAn deal with motion artifacts that are not uniform in the field of view?

CaImAn (Calcium Imaging Analysis) is a computational tool used for analyzing calcium imaging data, including correcting motion artifacts that can arise due to the movement of the sample during imaging. Non-uniform motion artifacts occur when different parts of the field of view (FOV) move differently, making the correction more complex than for uniform motion.

Describe the algorithm used to correct for elastic motion, and the steps it goes through.

To correct for non-uniform (elastic) motion artifacts, CalmAn uses an algorithm that involves the following steps:

1. *Partitioning the Field of View:* The FOV is divided into smaller overlapping patches. This allows the algorithm to independently estimate motion for different regions, accommodating local differences in movement.
2. *Rigid Motion Correction for Patches:* Initially, rigid motion correction is performed on each patch. This involves aligning the patches to a reference frame using translation and rotation adjustments.
3. *Elastic Motion Correction:* After the initial rigid correction, elastic motion correction is performed. This step uses non-rigid transformations to account for local deformations. Typically, this involves the use of techniques like optical flow, where the displacement vectors for each pixel or small region are estimated.
4. *Merging Corrections:* The local motion estimates from each patch are then merged to create a comprehensive motion correction map for the entire FOV.
5. *Applying the Correction:* The derived motion correction map is applied to the entire FOV to correct the recorded images, thus compensating for both rigid and non-rigid movements.

3. Constraints in the CMNF Algorithm

CNMF (Constrained Nonnegative Matrix Factorization) is an algorithm used to demix and denoise calcium imaging data by factorizing the observed data into spatial and temporal components that represent neuronal activity.

Constraints in CNMF:

- *Non-negativity:* Both spatial and temporal components are constrained to be non-negative. This reflects the physical reality that fluorescence and neural activity cannot be negative.
- *Spatial Smoothness:* The spatial components are often constrained to be smooth, reflecting the expectation that the fluorescence signal from a single neuron will be contiguous and not fragmented across the field of view.
- *Sparsity:* Temporal components are constrained to be sparse, meaning that neural activity (spiking) is expected to occur at discrete time points rather than continuously.

Reason for Constraints:

- *Biological Plausibility*: Non-negativity and smoothness ensure that the factors derived from the data make biological sense, as real neuronal signals and their spatial distribution are inherently non-negative and relatively smooth.
- *Noise Reduction*: Sparsity helps in reducing noise by focusing on significant, discrete events (spikes), thereby improving the robustness and accuracy of the inferred activity.
- *Improved Factorization*: These constraints guide the factorization process towards more meaningful and interpretable components, enhancing the quality of the decomposed signals.